N

Laboratoires Népendthès

AUTOTEST d'ovulation / Ovulation predictor self-test

PRINCIPE

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GB PRINCIPLE

- Upon release of an egg from an ovary, it passes through the fallopian tubes where it is ready to be fertilised. This is known as ovulation.
- Before ovulation, the body produces a large amount of luteinizing hormone (LH) which triggers release of a ripened egg. This phenomenon is known as the "LH surge" and usually takes place in the middle of the menstrual cycle. It is during this fertile period that a woman has the greatest chance of becoming pregnant. The test uses a combination of antibodies, including an LH monoclonal antibody which detects elevated levels of luteinizing hormone (LH).
- The Nep Ovulation Test detects the LH surge in the urine, indicating that ovulation will most likely take place in the next 24 to 36 hours. For pregnancy to occur, the egg must be fertilised by the sperm within 24 hours of its release. It is important to note that the LH surge and ovulation may sometimes be absent during a cycle.

PACK CONTENTS

- 5 test sticks • Leaflet • 5 sanitary gloves

PRECAUTIONS

Read the information provided in the leaflet carefully before doing the test.

- Do not use the test after the expiry date indicated on the foil pouch • Store in a dry place between 4-30°C. Do not freeze • Do not use the test if the foil pouch is open or damaged • Keep out of the reach of children • Test for in vitro diagnostic use only. Do not ingest • When you are ready to test, open the foil pouch • The test stick must only be used once.

INDICATIONS

The Nep Ovulation Test is for women finding it difficult to conceive and women taking ovulation inducers to confirm the effect of the prescription treatment.

The test is used to identify the period during which a woman is the most fertile.

BEFORE STARTING THE TEST

The kit contains 5 test sticks to enable you to detect ovulation on consecutive days. 1) Determination of the average length of your menstrual cycle. The total length of your menstrual cycle is the number of days from the first day of your last period (start of bleeding) to the last day before the next cycle. Use the table below to determine the day on which you should start the test.

Cycle length (days)	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38
Number of days before starting the test	6	6	7	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21

If your cycle lasts under 21 days or over 38 days, consult your general practitioner.

CALENDAR

Sun	Mon	Tu	We	Th	Fri	Sat
	1	2	3	4	5	6
7	8	9	10	11	12	13
14	15	16	17	18	19	20
21	22	23	24	25	26	27
28	29	30	31			

2) Circle the figure indicating the average length of your cycles on the table above. NOTE : If your cycles are irregular, you can refer to the shortest cycle. Morning urine should not be used to do the test.

3) Use the calendar below to determine the day on which you should start the test. For an optimal result, you should do the test at the same time each day. You should not drink in the two hours before the test. E.g. "My cycle usually lasts 28 days. My last period started on the 3rd. The table indicates 11 days to be counted as of the 3rd. The calendar below shows that I should start the test on the 13th day of the month."